User Manual | Bedienungsanleitung

visuReal[®] MASTER

HOYA
VISUREAL MASTER

visu^{GmbH} Solution

Dear customer,

thank you for purchasing the visuReal[®] MASTER video centration system!
You have made a good choice and we hope you enjoy using your new video centration system!

Please read the following information carefully! Ophthalmic knowledge is required for using this video centration system.

If you have questions, please contact your dealer. This is either the visuSolution GmbH or the respective HOYA subsidiary in your area! Our cooperation partner HOYA is authorised and well trained for supporting visuReal[®] MASTER.

Content

visuReal® MASTER

Content	04
Product Description	05
System Requirements	05
System Output	06
Mode of Operation	06
Technical Data	06
Scope of Delivery	06
Safety Instructions	07

First Steps

Assembly and Power Connection	08
Initial Startup	10

Operation

Measuring Process	17
--------------------------	-----------

Miscellaneous

Warranty	26
Disclaimer	27
Professional Disposal	27
Returning Goods	28

Product Description

The visuReal MASTER system is a novel video centration system, consisting of hardware and software components. It is a mirror system, containing a special camera system of six WQHD cameras. The extremely flat mirror case allows to be integrated into already existing or new shopfittings. The software is HTML-, e.g. browser based. The user can start the recording and measurement process and view the measurement results on a PC or notebook, tablet-PC or smartphone.

System Requirements

The visuReal MASTER is independent from any device and platform. That means it can be used on different end devices, on diverse operating systems and on common internet browsers.

Operating systems:

- iOS
- Android
- Windows
- Linux

Internet browser:

- Explorer
- Firefox
- Chrome
- Safari
- native Browser

End devices:

- Desktop computer
- Notebook
- Tablet
- Smartphone
- Large displays

System Output

visuReal MASTER automatically and precisely detects the centration data (PD, IPD, horizontal and vertical fitting data, pantoscopic angle, conreal vertex distance, head rotation, boxing, distance between the lenses, face form angle) by recording just one front-view image and without needing a measurement temple for attaching on the frame.

Mode of Operation

The customer stands in front of the mirror at an arm length's distance wearing the well adapted frame and looking into the mirror. The software is userfriendly and fast. The mirror provides you optimal support by evaluating the data automatically. All you have to do is to communicate with your customer, start the measuring process and receive the measurement data. The results can be checked and manually adapted, if necessary. There is no need for additional tools like a measurement temple or a side-view image.

Technical Data

- size: 64.5cm x 22cm x 3.5cm / 25.4inch x 8.6inch x 1.4inch (HxWxD)
- weight: 3kg / 6.6lbs
- camera: six-camera system, 8 MP WQHD
- power supply: 12V / 3 A

Scope of Delivery

- visuReal MASTER mirror unit
- wall mounting (wall mounting plate and screws)
- power supply (length 1.5m / 2m)
- user manual
- On/Off magnet switch

Safety Instructions

Remove the mirror always with both hands from the package to avoid possible damages.

The housing must never be opened. If doing so nonetheless, there is a risk of injury, e.g. through small damaged parts or mirror parts. The warranty will expire in case of an unauthorised opening.

Keep water or liquids away from the system. Do not expose the system to permanent jolts or vibrations.

Use only accessories which are part of the visuReal®MASTER-System (power supply and wall mounting), provided by the company visuSolution GmbH.

Do not use aggressive fluid solvents, chemicals or strongcleaning agents for cleaning the system. Use a soft cloth and commercial glass cleaner for the mirror surface.


The staff responsible for the operation and handling of the system must be informed about the mentioned safety instructions in connection with the system. Furthermore, they must be instructed to use the system correctly.

As the measurement system is supposed to be used with different end devices, please inform yourself about the respective safety instructions of your end device.

Assembly and Power Connection

Please note:

Please remove the device carefully from the packaging. Carefully remove the protective foil from the mirror side and watch out for dirt. When selecting the right place for the visuReal MASTER, please take care of a horizontally illuminated surrounding. Avoid backlight and direct spotlight on the person you want to measure.

Mirror Assembly:

Make sure, you align the wall mount correctly.

The mirror must be hanging **exactly vertically** to obtain **precise measurements**.

This means the mirror has to be 90 degrees to the floor and exactly vertically aligned.

Use the slots in the holder to align them with the wall.

Lay the power supply in such way that the cable is not under tension.

Power Connection and Automatic Startup:

Connect the visuReal MASTER with the 12V power supply:

The visuReal MASTER starts automatically.

When starting, the Eye Catcher status LED moves up.

As soon as the LED turns off (ca. 40 sec.) the system is ready to use. That means, it is ON when the light is off.

Assembly and Power Connection

Turn ON and OFF manually:

Turn OFF via magnet switch:

- The magnet switch sticks firmly on the bottom right side of the mirror
- Swipe the magnet to the left
- Then swipe it to the right
- The LED moves down for a few cm, stops for a moment, and finally moves down completely

Turn ON via magnet switch:

- Swipe the magnet to the left
- The LED-Licht moves up for a few cm, stops for a moment, moves up completely, then moves down to the middle and finally moves up.

How to find the right location in your shop:

- Make sure you have non-dazzling, horizontal light
- Avoid spotlight
- Make sure you have even illumination in the measuring area from all directions (from left and right side), so that one side of the face does not appear darker than the other side.
- When recording, take care, that there is nothing reflected in the visuReal Master mirror, which might be mistaken for the person that you want to have measured.
- Make sure, that there are no frame displays or many customers within the immediate measuring range of up to 80cm and behind it (up to 2m) (e.g. posters with faces, a waiting area or a cashier area with many people in the background may confuse the mirror system when shooting)
- Avoid backlight and reflections
- The visuReal Master can be hung on the wall between frame displays

Initial Startup

End Device Connection:

Connect your end device (tablet, laptop, PC) with the visuReal MASTER (initial start-up):

WiFi configuration

WLAN network: Master#
password: visuReal

Open the WLAN menu of your end device, search for the visuReal MASTER WLAN (a combination of the term MASTER and a serial number) and enter the password.

Starting the First Measurement:

NOTE (only for ninitia start-up + in case of changing the WLAN:

Please wait at least 3 minutes before starting the browser and entering the software URL.

Open the browser and enter the following URL:

URL for software

<http://192.168.77.30:8055/#/start>

Bookmark or link this open webpage to access it at any time.

You will always control the visuReal MASTER from this browser window. There you get all displays and results.

Troubleshooting:

If the system is not working properly, please check the following points:

- in case of initial start-up take care to wait 3 minutes before starting the browser and entering the URL
- The system is on (top LED is on). You can turn the mirror on and off by means of the provided magnet switch.
- WLAN must be connected to MASTER#, connection must have sufficient signal strength.
- Check the URL for typing errors.


Initial Startup

Connect visuReal MASTER into an Existing Network (optional):

Open your browser and enter the URL for visuReal MASTER.

Click on the hamburger menu and select WiFi settings.

Select your respective WLAN and enter the respective password.

Network settings for internet access and software updates

Softwareupdate via VPN
only when internet access is available

Available WiFi Networks
Please select a WiFi providing internet access

SSID

WiFi Key

key

back

save

Initial Startup

AutoUpdater:

Getting updates becomes easy:


Thanks to the AutoUpdater-feature you will receive an automatic update notice in form of a red one (1) in the hamburger menu as soon as there is an update available. If you click on it, a security prompt appears and then the update is carried out.

NOTE If you want to use the AutoUpdater feature, you have to make the following settings at first:.

Go to the settings of your end device that you use for visuReal Master, e.g. your iPad, and go to „Network settings for Internet Access and Software Updates“
 Activate the slider

- Updates and Service via VPN
- Automatic information about available updates

Once these settings have been made, the system checks the central update server each time it boots up, i.e. App start, to see whether an update is available.

Initial Startup

Automatic WLAN- channel optimisation of the internal 2.4GHz-Master hotspots:


WLAN issues are an old hat:

In case of an inadequate running WLAN connection between the end device and the visuReal Master (the mirror) there is a solution now:

For this, please go to the hamburger menu in your visuReal Master software and go to „Network settings for Internet Access and Software Updates“
Activate the setting „Automatic channel optimization“.

This switches the WLAN to a channel, that is least affected by the number of other hotspots. This can be done automatically each time the program is started (switch on „Regularly at system start“) or once by pressing the [Once only]“ button.


Please note: Here, it is necessary that the WLAN settings of the end device are set to „connect automatically“, otherwise the WLAN connection to the visuReal Master must be re-established manually.

Initial Startup

Export Settings for Measurement Data Export and E-Mail Settings:

Go to the hamburger icon and select the export settings. There, you enter the data required for exporting the measurement data lateron:


Data entry for the data export to HOYA:

Please fill the fields customer number, login and password of your Hoya-login data.

Data e-mail:

For easily transferring measurement data by e-mail, enter the e-mail address to which you always want to send your measurement data into the field at the bottom. This may be the e-mail account of your device on which you are using visuReal Master or any other.

Initial Startup


Software Settings:

Before starting with the first measurement, you can configure the software. For this, click or tap on settings icon for above in menu bar.

Automatic exposure time:

The checkbox for „Auto-Exposure“ is set by default.

Contrast:

The default value is 0.6

The contrast can be adjusted infinitely by moving the slider, if it is not satisfying.

Exposure time:

The default value is 0.2

In case of too bright or too dark images the exposure time can be extended or shortened by moving the slider.

Eye-Catcher:

you can activate or deactivate the Eye-Catcher by default by using the slider (default setting ON).

Show preview images:

Here you enable or disable the display of the measurement images for the measurements listed under History of the last 48 hours (default setting ON)


Initial Startup:

Other Settings:

History:

In the hamburger menu you will find the option History. All measurements are stored for 48 hours and can be viewed under History and edited again. The edited measurement is saved again under the same order number for 48 hours.

After 48 hours, the measurements are automatically deleted. You can also delete them manually.


Language:

It is possible to switch to another language.

Click or tap the flag icon next to the hamburger menu. an overview with all available languages appears.

The default is always English. Click the language you want to use.

Measuring Process

Important: Adjust the selected frame correctly to your customer before starting the measurement! To avoid unintended reflections, please remove the demo lenses, if possible. Take care of a natural head and body posture.

Record Image

Your customer stands in front of the mirror at arm's distance (70 to 80 cm/ 27 to 31 inch) wearing the well adapted new frame. There is no need to put a measurement temple on the frame. The customer just looks in the mirror at his nasal root or at the red LED of the Eye Catcher.

Please note about the Eye Catcher:


Measuring with Eye Catcher:

Customers, who cannot see themselves in the mirror without spectacle lenses, need the Eye Catcher's red LED to for focusing (default setting, recommended).

Measuring without Eye Catcher:

Customers, who can see themselves in the mirror, do not need the Eye Catcher. They look in the mirror focusing their nasal root.

Each measurement is automatically provided with a time stamp. Optionally, you can enter an order ID for a better allocation.

Please press **[Start]**.


The mirror detects the person and detects automatically, which 2 of the 6 cameras are supposed to be used for the recording.

As soon as the **[Shot]-button** is active, press or tap the button in order to trigger the recording. A front-view recording is done, processing two almost identic images.

Check image


After a few seconds, you will see the two combined stereo images, that allow a 3D evaluation. Check the centration marks. There are ten (10) cross marks for the frame and two (2) iris and pupil cursors for the eyes in the top image. In the related bottom image there are eight (8) cross marks for the frame and two (2) iris and pupil cursors. If these markings are not correctly recognized, each marking can be tapped, marked and newly positioned.

If everything is correct, please click or tap on **[results]** to access the result page.


Explanation markings

Basically the following applies: The crosses for the frame have to be set on the most outside point on the inner frame edge respectively at the top, at the bottom, at the nasal and temporal side of the frame edge. The circle cursors for the eyes have to edge the iris respectively the pupil. The crosses marked in red in the illustration, have to be placed as far as possible sideways to the temporal or nasal side just before the frame is curving down. The temporal and nasal crosses marked in green you place nasally resp. temporally at the extreme point of the frame on the inner frame edge. All markings are blue in the software before they are manually corrected. The red and green circles are for explanation only in this document.


Always check the centration marks. If needed, adapt the marks manually. Click or tap the crosses you want to correct. Once selected, the blue crosses turn red. Then press the button **[correction]**.

Correct markings


Now, for each selected marking, there will appear a correction program step. Now you can move the the red little flags with your mouse or fingertap (depending on your end device). If you have selected several crosses, they will be displayed subsequently in single images.

First the same crosses, then the cursors for the eyes. You pass to next correction image with **[next]**. If everything is done, you can complete the process and go to **[recalculation]**.

Please note the magnifying glass view

For better usability, there is the magnifying glass displayed. In order to be able to place the measuring points easier and better with the mouse or the finger, a magnifying glass view appears when „touching“ the measuring points. If the measuring point is moved again after a manual correction, the loupe circle of the starting position remains visible. Thus, you know where the measurement point originally was (see image 5, magnifying glass view)

Please note the cursor connection:


The markings of the two images are connected to each other. When moving the cursor in one image, the same cursor will move in the other one. This does not apply for the green nasal and temporal crosses illustrated on page 16. For them, there is only the top image displayed in the correction mode.

Measuring Results

After clicking on [results] or [recalculation], you will pass on to the final program page showing all results at a glance:

Important! When terminating the measurement process, the results are stored for 48 hours. Then, those measurements are automatically deleted (See page 14 Other Settings - History).

Additionally, to the text view there is a graphical illustration showing you the results. This makes it easy for you to assign the values quickly, to perceive them at first sight and to assess whether corrections are necessary or not.

You will see the following data at a glance.

Organisation Data:

- order ID (if entered)
- time

Centration Values:

- IPD
- monocular PD
- fitting height (y to frame)

Head Posture:

- back vertex distance BVD
- inclination

Frame Dimensions

- boxing (A+B size)
- Distance betw. lenses DbL
- Frame type (must be selected)
- Face form angle


Lens:

- Lens type (must be selected)

Select Frame Type

On the result page, you have to select the frame type, i.e. the frame material. Please select according to the frame you have chosen:

- Metal
- Plastic
- Rimless
- Nylor
- Optyl
- SPX

The frame type selection is used for calculating the groove depth of the eyewire, which is integrated into the centration data and transferred to Hoyailog during the export.

Important! Please select the frame type. This is obligatory for receiving correct results! At the end, the measurement results are adapted respectively.


Select Lens Type

On the result page, you have to select the lens type and the respective centration method.

Please select between:

- progressive
- Single-Vision Pupil Center
- Single-Vision Center of Rotation

The measuring results are provided according to the selected lens type. In case of single vision lenses according to the center of rotation, there will be different results for the fitting heights $y(R)$ and $y(L)$. In case of progressive lenses, the largest possible progression zone lengths will be determined and displayed.

Calculation of the head rotation correction

You have the possibility to compensate the measured head rotation, if you think that it is needed. The head rotation is displayed in degrees. If you want to correct the measured value to 0 (zero) degree, you have to activate the slider [rotation compensation] on the measurement result page. The first value shows the head rotation in degree (here in the picture 4.5). The second value (here in the picture 2.2) indicates by how many millimeters the monocular PD has changed, having activated the [rotation compensation] to zero degree. In the following you can see two program screenshots showing the correction to zero, i.e. "rotation compensation" OFF and ON.:

Rotation rotation OFF

Centration Values		
interpupillary distance:		64.5
x	R 31.9	L 32.6
y (to box)	R 18.3	L 17.3
Head Posture		
back-vertex-distance:		12.1
inclination:		13.5
rotation (° mm)		4.5 0.0
<input type="checkbox"/> rotation compensation		

Rotation rotation ON

Centration Values		
interpupillary distance:		64.5
x	R 34.2	L 30.3
y (to box)	R 18.4	L 17.4
Head Posture		
back-vertex-distance:		12.1
inclination:		13.5
rotation (° mm)		4.5 2.2
<input checked="" type="checkbox"/> rotation compensation		

Manage Measurement Values

On the result page, you have the following options to further process your measuring data.

Button [edit]:

Here, you can once again go back to the correction page and edit the measurement by tapping the button [edit].

Button [data e-mail]:

Here you can save the results on your end device. Tap the button [data e-mail] and send the measuring results to an e-mail address. For this, you have to enter the desired e-mail in the export settings (see page 12 export settings).


Button [export data]:

Or you can export the data to the HOYA order system. When tapping the button [export data], the following data are transferred to HOYA. The calculated values already contain the information from the selected lens and frame type.

- IPD
- monocular PD
- fitting heights
- distance between lenses
- Face form angle
- back vertex distance BVD
- inclination
- A size
- B size
- order iD as reference

Button [finished]:

Here, you finish the measurement procedure and go back to the start screen.

„inSession“ Feature for the data transfer to HoyaiLog

You can now add the visuReal Master values to an already existing data record (inbox). Thus, you can manage your data even better.

When touching the button data export in order to export to Hoya, there is now an additional switch (see marking in image 6). By means of this switch, you confirm to add the measurement results to an already existing dataset of the inbox. To do so, turn the switch ON.

Then, the inbox list appears. From there you select the required dataset to which you want to add the visu-Real Master measurement results. Then click on YES.

After confirming and successfully uploading the data, the following message appears “saved successfully”

Submit order

Are you sure, you would like to export your measurement data?

Add measurement to an existing consultation

No
Yes

Submit order

Are you sure, you would like to export your measurement data?

Datum	Programm	Kommission
2019-11-22 15:08:15	master	Test_Hendrik_20191122_1507
2019-11-22 13:05:33	HoyaiLog	Hendik
2019-11-22 13:02:49	MyStyle iDentifier	Hans Joachim Ollendorf
2019-11-22 13:02:41	HoyaiLog	Kunststoff Breitbügel
2019-11-22 13:02:38	HoyaiLog	toth mátyás

Add measurement to an existing consultation

No
Yes

Warranty

The warranty period granted for the visuReal MASTER video centration system is 12 months. The period of warranty shall begin on the date of delivery. Services under warranty will be provided by the company visuSolution GmbH or their sales partners. In case that you notice a malfunction of the product, immediately give a notice to the company visuSolution GmbH (see address below) or to their sales partners.

The warranty is only valid in countries where the visuReal MASTER video centration system is sold. There is no on-site warranty, that means the user will pay possible transport costs.

The company visuSolution GmbH shall not assume liability for transport damage if the user did not complain about visible package damage upon receipt.

If the product becomes defective within the warranty period, it will be repaired or replaced free of charge, provided that it has been installed and used properly and unless it has been modified or repaired by the user. The warranty will not be valid if the product has any defect caused by misuse, incorrect handling, storage or modification.

The warranty will also expire if the mirror housing has been opened by the user. Additionally, the warranty will expire automatically if the type plate/serial number were detached from the hardware components.

Services provided under warranty do not extend or restart the warranty period. Replaced defective parts or apparatus will become the property of the company visuSolution GmbH.

If the cause of damage has been improper handling, misuse, or non-observance of the safety instructions or the operating instructions, the user shall be obliged to cover the costs for the warranty service.

The product visuReal Master leaves our premises after a hundred percent final end inspection. We have high quality standards and define them much higher than our suppliers. Minimal inclusions or scratches on the mirrored surface are due to machine production. These in no way affect the function or the measurement results and are therefore excluded from a complaint. As a reference condition in this regard, a distance of 60cm to a vertically oriented visuReal Master is determined. Please note the chapter Assembly under this document and the assembly sketch included as a hardcopy in the delivery.

Disclaimer

The company visuSolution GmbH shall not be liable for the compensation of any special, indirect, incidental or consequential damages caused by the use or the improper use of the visuReal MASTER system. This exclusion of liability includes, but is not limited to loss of profits, loss of investments, lost revenues, sales or data (whether direct or indirect) or commercial loss for breach of any expressed or implied guarantee on the visuReal MASTER system, even if the company visuSolution GmbH has been advised of the possibility of such damages. To the extent permitted by applicable law, this shall also include damage to property and personal injury.

Except to the extent prohibited by applicable law, any implied warranty or condition of merchantability or fitness for a particular purpose (see product description) of the visuReal MASTER system is limited in duration to the duration of the warranty period.

Some states do not allow limitations on how long an implied warranty lasts, or do not allow the exclusion or limitation of liability for incidental or consequential damages, so the above limitations may not apply to these states. This warranty gives you specific legal rights. You may also have other rights that vary by country or other jurisdiction.

The company visuSolution GmbH shall not assume liability for damages arising directly or indirectly from the incorrect use of the visuReal MASTER system due to translation errors in the operating instructions (user manual). The company visuSolution GmbH shall not be liable for damages directly or indirectly caused by incorrect descriptions, omissions, changes without prior notice, different interpretations or differences between the operation of the system and the descriptions in the operating instructions.

Professional Disposal

This product is subject to European Directive 2002/96/EC. This icon means, that this product must not be disposed of with general household garbage. It shall be returned for recycling. In case of questions, please contact your seller.

Returning Goods

In case you have purchased your system from HOYA, please contact your respective HOYA subsidiary in the first instance in the event of a complaint / return consignment.

In case you have purchased your system directly from visuSolution GmbH, please contact the visuSolution GmbH.

To ensure a quick completion of your request, it is necessary to follow the steps below:
In case of a complaint/ return consignment please request an RMA number at:

**telephone: +49 39361 967-17 or
e-mail: rma@visusolution.com**

When returning your consignment please use the enclosed RMA sheet or attach the RMA number well visible outside on the parcel.

Return consignments without RMA number will not be accepted by us and will be sent back to you!

Replacement deliveries or credits notes will be effected promptly after verification:

1. The return consignment has to be carried out carriage-paid.
2. Please fill out and attach the RMA sheet completely including a detailed error description.
3. Please always use the following delivery address for your return:

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4. Please always use shipping package for packing the products. Do not label or write on the original packaging.
5. Please return goods within 10 days after requesting the RMA number. Goods that arrive later will be rejected.

The return consignment cannot be accepted, if one of the previous points were not considered. This return consignment is subject to our general terms and conditions which can be found on our website.

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www.visuSolution.com

designed and developed by  **Ollendorf**

|visuReal[®] is a registered trademark of visuSolution GmbH

copyright March, 2019 visuSolution GmbH, edition 9, software version 1.3
Technical changes and errors reserved

photography copyright: HOYA Holdings for visuReal MASTER eyecare professional and product shots

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Copyright 05. März 2020 visuSolution GmbH, 9. Auflage, Softwareversion 1.3